ROOYAGEN

Platelet Rich Plasma Gel Kit

Platelet Rich Plasma Gel Kit

CHRONIC WOUND HEALING

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
REFERENCE

- Mustoe T. "Dermal ulcer healing: Advances in understanding". Tissue repair and ulcer/wound healing: molecular mechanisms, therapeutic targets and future directions. Paris, France: EUROCONFERENCES, March 17-18, 2005
- Snyder, Robert J. "Treatment of nonhealing ulcers with allografts". Clinics in Dermatology (2005). 23 (4): 388-95.
- Augustin, M.; Maier, K.. "Psychosomatic Aspects of Chronic Wounds". Dermatology and Psychosomatics(2003);4:5.
- Velander, Patrik E.; Theopold, Christoph; Gheerardyn, Raphael; Bleiziffer, Oliver; Yao, Feng; Eriksson, Elof. "Autologous cultured keratinocytes suspensions accelerate re-epithelialization in the diabetic pig". Journal of the American College of Surgeons (2004). 199 (3): 58.
- Supp, Dorothy M.; Boyce, Steven T. "Engineered skin substitutes: Practices and potentials". Clinics in Dermatology(2005). 23 (4): 403-12.
- Centers for Disease Control Diabetes data and trends. Available at: www.cdc.gov/diabetes/statistics/complications_national.htm cessed January 23, 2012.
- Kumada Y, Ooura T. February 10 as "Foot Care Day". Office of Diabetes Net. Available at: www.dm-net.co.jp/calendar/2012/016573.php. Accessed January 16, 2012
- CMS letter for non healing chronic wound by platelet-gel, CAG-00190R3, feb, 2012
- Alsousou J1, Ali A, Willett K, Harrison P. The role of platelet-rich plasma in tissue regeneration. Platelets. 2013;24(3):173-82
- Driver VR1, Hanft J, Fylling CP, Beriou JM; Autologel Diabetic Foot Ulcer Study Group. A prospective, randomized, controlled trial of autologous platelet-rich plasma gel for the treatment of diabetic foot ulcers. Ostomy Wound Manage. 2006 Jun;52(6):68-70, 72, 74 passim
- de Leon JM, Driver VR, Fylling CP, Carter MJ, Anderson C, Wilson J, Dougherty RM, Fuston D, Triglia D, Valenski V, Rappl LM The Clinical Relevance of Treating Chronic Wounds with an Enhanced Near-Physiological Concentration of Platelet-Rich Plasma Gel. Advances in skin & wound care 2011 Aug 24 8 357-68
- Wagner Grading System for Diabetic Foot Infections, Original Text by Clifford R. Wheelless, III, MD. Last updated by Data Trace Staff, Dec- 2012.
- Bennett NT1, Schultz GS. Growth factors and wound healing: Part II. Role in normal and chronic wound healing. Am J Surg. 1993 Jul;166(1):74-81.
- Sheehan P1, Jones P, Caselli A, Giurini JM, Veves A. Percent change in wound area of diabetic foot ulcers over a 4-week period is a robust predictor of complete healing in a 12-week prospective trial. Diabetes Care. 2003 Jun;26(6):1879-82.
- FDA cite: 21CFR640.25, part C, platelet preparation, WWW.FDA.gov
- Dougherty EJ, An Evidence-Based Model Comparing the Cost-effectiveness of Platelet-Rich Plasma Gel to Alternative Therapies for Patients with Nonhealing Diabetic Foot Ulcers, Advances in Skin & Wound Care: 2008 - Volume 21 - Issue 12 - page 568-575
- Wu SC, Driver VR, Wrobel JS, Armstrong DG. Foot ulcers in the diabetic patient, prevention and treatment. Vasc Health Risk Manage. 2007;3(1):65-76
- Senneville, Eric MD 1; Lombart, Audrey MD 1; Beltrand, Eric MD 2; Valette, Michel MD 1; Legout, Laurence MD 1; Cazaubiel, Marie MD 1; Yazdanpanah, Yazdan MD, PHD 1; Fontaine, Pierre MD, PHD. Outcome of Diabetic Foot Osteomyelitis Treated Nonsurgically: A retrospective cohort study. Diabetes Care. 31(4):637-642, April 2008.
- Lipsky, B. A. 1, 2 Empirical therapy for diabetic foot infections: are there clinical clues to guide antibiotic selection? Clinical Microbiology & Infection. 13(4):351-353, April 2007
- International Cellular Medicine Society - 2011, section VIII Platelet Rich Plasma (PRP) Guidelines www.cellmedsociety.org
- Treatment of osteoarthritis of the knee, evidence based guideline 2 ND edition, AAOS board of directors, may 18/2013
- FDA 510(k) clearance, BK060007. AutoloGel System. Rockville, MD: Cytomedix, Inc; September 20, 2007
- Anjayani S, Wirohadidjojo YW, Adam AM, Suwandi D, Seweng A, Amiruddin MDInt J Dermatol. Sensory improvement of leprosy peripheral neuropathy in patients treated with perineural injection of platelet-rich plasma. 2014 Jan;53(1):109-13.
- Sakata J, Sasaki S, Handa K, Uchino T, Sasaki T, Higashita R, Tsuno N, Hiyoshi T, Morimoto S, Rinoie C, Saito N. A retrospective, longitudinal study to evaluate healing lower extremity wounds in patients with diabetes mellitus and ischemia using standard protocols of care and platelet-rich plasma gel in a Japanese wound care program?. Ostomy Wound Manage. 2012 Apr;58(4):36-49.
- Schwartz principles of surgery text book, 10th edition, 2014
- Mohammadi, S. Nasiri, S. Mohammadi, M.H, Malek Mohammadi, A, Nikbakht, M. Zahed Panah, M, Safar, H, Mostafaei, S, Norooznezhad, A.H, Soroosh, A.R, Alimoghaddam, K. Ghavamzadeh, A. Evaluation of platelet-rich plasma gel potential in acceleration of wound healing duration in patients underwent pilonidal sinus surgery: A randomized controlled parallel clinical trial. Transfusion and Apheresis Science. Volume 56, Issue 2, April 2017, Pages 226-232
- Mohammadi, M.H., Molavi, B., Mohammadi, S, Nikbakht, M., Mohammadi, A.M, Mostafaei, S, Norooznezhad, A.H., Ghorbani Abdegah, A., Ghavamzadeh, A. Evaluation of wound healing in diabetic foot ulcer using platelet-rich plasma gel: A single-arm clinical trial. Transfusion and Apheresis Science. Volume 56, Issue 2, April 2017, Pages 160-164.

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When pulling the upper plasma in the first stage, because of the Negative pressure inside the tube, which leads to the plasma return from the syringe to the tube, insert a sterile needle beside the IV Cannula into the tube.

3- Then pick three B tubes in the kit and through the adaptor without any pressure on the piston of the syringe, transfer the blood into them and fill them up to the upper edge label on the tube (to fill the one cc at the end of the tube, avoid pressing the piston of the syringe and the blood should be entered into the tube through inverse pressure, in a way that the piston is pulled back and released to enter the blood into the tube by evacuation).

4- Put 3 tubes B in balance into the centrifuge for 10 minutes adjust its turning in accord with the following table.

5- After centrifuging, through peripheral IV cannula 16G and 5ml syringe inside the kit which are in stage II, transfer the upper plasma and the buffy coat into tube C.

6- In total the plasma in upper parts of 3 tubes B will be transferred to 3 tubes C (to pull some RBC along with the buffy coat is not a problem).

7- Then, add 2.5ml of the activation solution to each tube and wait 20-30 minutes in room temperature for the gel formation (to accelerate the formation to gel, the incubator heat in 37°C or hand heat can be used).

8- In each tube one roll of 10 ml volume will be formed. By each roll a wound having 6-8 Cm² can be covered. In the end, by one kit containing 3 rolls gel an area about 3 times as much the above size can be covered.

9- Then, put the gel on the wound and dress it with the Vaseline gauze in the kit without any pressure.

Use the safety box to remove materials, needles, syringes, plasma and components.

Brand	Model	Status	Country	10 min (PRM)
Rooyagen	-	Swing	Iran	1600
Hettich	EBA 20	Fix	Germany	1200
Hettich	EBA 270	Swing	Germany	2700
Hettich	Rotofix 32	Swing	Germany	1600
Pars azma	-	Swing	Iran	1500
Jenus	TDL-4	Fix	China	1000
Sehand	-	Swing	Iran	1600
Labtron	-	Swing	Iran	1600
Selecta La	-	Fix	S. Korea	900
Selecta La	-	Swing	S.Korea	1300

Side Effects:

In order prepare platelet gel; the individuals' own blood will be used. Thus, it has no side effect necessarily unless the tubes, containers of inferior quality, laboratory solutions and unauthorized kits produced from industrial plastic substances (not medical) have been used.

Executive indices and expectable results:

81.3% of the wounds lower the 7 centimeters in average time between 6.2 weeks in a random clinical trial (RCT) for diabetic wounds of foot were healed completely.⁹⁻¹⁰

In a large clinical trial, an average 15-day improvement in subjects treated with PRP gel was significant compared to a mean of 337 days in the control group without this treatment. And 96.5% of all ulcers with an average of 2.2 weeks and an average of 2.8 times of dressing in 285 patients that were observed in 39 centers.¹¹

In a longitudinal prospective clinical study, with regard to the increase of depth (3.5x), area (8.8x), and volume (16.3x), the wounds, after 75 days of standard treatment, with the beginning of platelet gel treatment with the average of 6.1 times of dressing had 83% perfect healing and the advancement of patients to later phases and amputation of lower organs were prevented.

In one of these studies, Wound healing after pilonidal sinus surgery with and without Rooyagen PRP gel was evaluated. According to the results, patients treated with PRP gel went through a significantly faster healing process (8.69 ± 1.18 in controls and 4.78 ± 0.87 weeks in PRP gel treated ones with the P-value = 0.03) and returned to their routine activities (3.3 ± 0.64 for the treatment of arm and 6.5 ± 1.03 weeks for controls with the P-value = 0.00) while experiencing less pain (P-value = 0.00) and shorter anti-biotic consumption duration (P-value = 0.00). Conclusion Considering the results, authors of this study suggest PRP gel treatment for post operation wound dressing of pilonidal sinus disease with healing by secondary intention²⁶. And in another study investigators evaluate the effectiveness of using autologous Rooyagen PRP-Gel for treatment of diabetic foot ulcer (DFU) during the first 4 weeks of the treatment, were shown PRP could be considered as a candidate treatment for non-healing DFUs as it may prevent future complications such as amputation or death in this pathological phenomenon.²⁷

- Sickle Cell Anemia Wounds
- The wounds resulting from spinal cord injuries
- The clinical infection wounds

Contra indication:

- Patients under treatment by chemotherapy drugs
- The wounds caused by malignancy
- The Active clinical infectious wounds
- Non-consumption of NSAID drugs and chelating agents for calcium one week prior to treatment
- Lack of disorder in primary homeostasis (platelet syndrome) and secondary homeostasis (coagulation factors)
- Non-consumption of caffeine, garlic, ginger, omega 3, and omega 6, one week prior to treatment
- Systemic disease such as diabetes with weak control of HbA_{1c}>8 % uncompensated cardiovascular diseases, the weakness of immune system and taking drugs suppressing immune system
- Taking anticoagulants and drugs disturbing platelets performance such as Ticlid and Plavix
- The count of platelets less than $10^5/\text{mm}^3$ (100000 in a microliter) and the anemia Hb<10 gr/dl
- The patients infected by topical abscess or arrhythmia and sensitivity was excluded of treatment.

Explanation and warnings:

- The blood drawing should be done by experienced individuals without any stress, in one time period (about 3 minutes).
- It is preferred to pour the blood through the walls into the tubes (without any extra movement which can lead to damage to platelets).
- Make sure of the balance of centrifuge and the area between blood and plasma should be sharp in the end. If this area is not seen, the balance of the table under the centrifuge and its tightness and, finally its problem, should be assessed. After centrifuge the concentration of platelet should be 2 ± 0.2 times as much of that of the basic platelet which is sufficient for the production of the gel.
- Do not change the components and accessories associated with the different Batch numbers.
- The contents of this kit are for one use and a patient.

Kit contents:

Stage I

Syringe: The 30 ml Volume Syringe For Sampling

Sampling Needle: Butterfly Needle (21G) for sampling (1X)

Tube (A): Tube a containing 5ml anticoagulation and blood cells separator (1X)

Adaptor: Adaptor for transferring blood from the syringe to tube B (1X)

Tube (B): 10ml tubes with non- adhesive walls to platelets (prp separation tubes) (3X)

Method: The method of work will be in according to ifu (brochure of the kit)

Stage II

5ml Syringe: 5 ml syringe to carry samples (2X)

Needle 21 G : For transfer (1X)

IV cannula: (16 G) for transferring plasma (1X)

Tube (C): With 10 ml volume containing 1 ml of gel stabilizer (3X)

Activation vial: Containing 10ml vial for gel activation reagent (1X)

Sterile Vaseline gauze: For dressing after topical gel application (2X)

The material and devices not found in the kit:

Centrifuge for 16×100 tube; swinging arms preferred
Safety Box
Disposable Gloves
Tube Cases (Rack)
Alcohol 70% Spray
Tourniquet
Forceps

Sustainability and maintenance conditions:

The kit and its components can be used up to the expiry date written on each of the components provided that they have been packaged and their package has not been opened and damaged.

The storage temperature of the kit is preferably room temperature (RT), but not kept below 2 ° C or more than 30 ° C

Procedure:

PRP gel kit will be prepared through making use of 27cc of the Peripheral blood and within 10 minutes the PRP with the dosage in mind will be prepared as follows:

1- At first, the stage I of the kit should be opened and "3cc of solution A" with 30cc syringe should be pulled and through putting scalp vein needle, 27cc of the blood should be drawing in order to the volume of the blood and the solution inside the syringe reach 30cc.

2- Connect the adaptor of blood transfer to the syringe and upside-down the connected syringe four times slowly.

pulsed lavage will be sufficient. In most cases, the cleansing of wounds with physiologic serum and keeping the wound wet with this serum will be helpful, for healing of the wounds lower than 3 cm size. But, in case of flesh Tissue wounds which have been classified as non-healing wounds and the other deep and complex wounds (as explained in gradation part), especially in diabetic patients infected by the wounds resulted from advancing diabetes specially in lower organs (foot) (non-ischemic cases) which have been categorized as non-healing wounds, the dressing of wound through making use of the platelet gel is a helpful method.¹⁷

In case of the presence of abscess, it should be evacuated and the wound should be debridement. In case of deep wounds in which tendons and joint capsules are involved, the debridement of the wound and the removal of necrotic tissues are necessary and, after that, the plastic surgery may be needed. Non infected cellulite wounds resulting from deformation of foot or due to neuropathy do not need antibiotics after debridement and it is possible to treat them by occlusive cast.¹⁸

In cases where the wound is deep and the abscess is formed and the cellulite is present, or there is infected arthritis and osteomyelitis, prescription of antibiotics is a must. Before antibiotics have been prescribed the aerobiosis and anaerobiosis culture and antibiogram of the wound should be performed. The sample should not be taken from the surface of the wound, rather. It should be obtained from the depth of the wound until the result of antibiogram is known, a Broad spectrum antibiotics which affects both positive gram and negative gram bacteria should be prescribed.¹⁹

As the majority of the microbes resulting in infection are positive gram cocci, that is, Staphylococci and Streptococci, an antibiotic to which the microbes are sensitive should be prescribed. It is better to prescribe antibiotics such as Cephalexin and Carbapenems for the treatment of such patients. Of course, some of the authors have reported good results for the prescription of oral Clindamycin in 750 milligrams per a day. This antibiotic has proved to affective on positive gram cocci, negative gram strains, and 3/4 of anaerobiosis microbes. The period of prescription depends on the severity and the depth of the damage. In cases where the infection has not penetrated to the bones, two weeks antibiotic prescription will be sufficient which, depending on the conditions of the patient, can be used by eating or drinking. In some cases in which the bones have been infected a longer period of treatment is needed.²⁵

In a study performed through making use of PRP-Gel ROOYAGEN, the patients' wounds were bandaged topically by platelet rich plasma gel and dressed by Vaseline gauze available in the kit. After 24 to 48 hours

(preferably 36 hours) the dressing of the gel will be removed and, after washing with normal saline the standard cares will be done up to the next referring (next week). Although the signs of healing can be seen from the second week and the average period need for the perfect treatment lasts from 6 to 8 weeks, most of resources do not advise more than 12 weeks dressing of the wound. The evidences of clinical response as the capability of generating positive results, intensity of response, the time of response and the number of medications have been taken into consideration to reach a positive response.

Dressing with platelet gel will have a significant influence on the quality of life's of the patients through the reduction of the costs and expenses which prevent from amputation and shorten the healing period has always been taken into account for the cases resistant to popular treatments. Platelet rich plasma gel contains cytokines, growth factors, and chemokine and fibrin membranes from the patients' blood. The mechanism of action of PRP gel in cellular and molecular induction of the wounds' repairs response similar to that seen in platelet activity.

The most important challenge facing PRP is separating rich platelets from the blood with controlled concentration for different indications. In this preparation the requirements of FDA (No. 21CFR640.25) should be admitted. In addition to platelet swirling, the platelet count by CBC test should be confirmed by the laboratory (for at least 1% of the referrers).

Where PRP gel is used the most and where not used

Today, several special PRP kits have received the confirmation of FDA (American food and drugs) and the CE of European Union and there are more than 10000 articles indexed in national center for biotechnology information site (NCBI) on how to use PRP gel in treating non-healing chronic wounds such as diabetic wounds, skeletal-muscular, surgical usages, hair and skin beauty, dentistry and the like²⁰. In addition to these, this technique has been able to receive the license of CAG-00190R3 for the treatment of chronic wounds⁸, the confirmation of FDA for dressing chronic wounds, and the agreement of American Association of Orthopedic Surgeons (AAOS) in the report in 2013.²¹⁻²²

Indications:

- Diabetic wound of foot
- Pressure wounds
- Vascular wounds
- Traumas
- Dehiscence wounds

of Medical Care Service in 2012, the results of 25 years research have been published in 100 pages which support the effectiveness of this claim.

⁸ In addition to the abundant studies in this field, more than 21 random clinical trials have also been published. In these studies, in addition to the references used in the current study, three systemic studies or clinical trials which have considered all of the aspects of comprehensive study and the information related to the perfect or relative healing of the wound, time length and the procedures of the healing, the rate and acceleration of the repair, and the reduction of the wound size have been taken into consideration in these all-inclusive studies.⁹⁻¹¹

There is various classifications for diabetic wounds, the most famous and simplest of them which are still being used in the Wagner scoring system where the diabetic ulcers of the foot have been graded in 5 phases:¹²

Group 0- In this phase the skin of the foot is intact and unchanged but a deformation exists in foot bones due to diabetic neuropathy and the increase of pressure on some parts of the foot make it susceptible to ulcers.

Group 1- Topical and flesh wounds

Group 2- A deep wound which, in addition to the skin and the hypoderm, involves the tendon of the joint capsule

Group 3- The affection the bones, gangrene of the toes and the front part of the foot

Group 4- Gangrene of the toes and the front part of the foot

Group 5- Gangrene of all parts of the foot

In another prevalent classification the diabetic infections of foot have been divided into two groups:

1-Surface infections: include nail bed infections, nail drooping, infections between the fingers and foot ulcers. If the flesh wound infections are not treated in due time, they will lead to deep infections.

2-Deep infections: These include the deep cellulite of the foot, the abscess of the sole and back of foot. In addition to the skin and the hypoderm, the deep infections involve tendons, joint capsule and the bone of the foot.

Topical and surface infections do not usually have any general symptom, but, the deep cellulitis are usually accompanied with ague and the inappropriate general conditions. Due to the long-term healing and in some cases the lack of improvement of secondary wounds to diabetes, which leads to amputation, especially in the wounds of the lower extremity, and in the event of failure in the treatment of surgery, amputation.

Types of chronic wounds with regard to recovering¹³

Healing wounds:

- The growth factors derived from the platelets, cytokines, and chemokine are present in time and place and in a balanced manner
- There are sufficient perfusion in order to carry micronutrients to the wound and sufficient angiogenesis in location of the wound
- The processes of reconstruction, structure of the wound and the tissue matrix have a natural progress.

Non-Healing Wounds

- Imbalance in the growth factors derived from the platelets, cytokines, and chemokine disturb the process of healing.
- The structure of inflammatory cytokines leads to the creation tissue destruction cycle.
- The lack of an appropriate blood perfusion and its support limits the feeding of the wound
- There is not any progress in reconstruction, structure of the wound and tissue matrix.

Complex Wounds

- Tubular wounds:** One of the problems of these kinds of wounds is the diagnosis of them within the tissue. As most of the damages happen under the skin and not visible, the healing of this kind of wounds last longer than the normal ones.
- Sinusoidal wounds:** These are wounds to which the access is difficult.
- The wounds with visible organs** including blood vessels, nerves, leprous wounds, tendons, ligaments, and/or bones²³
- Bedsore:** These are the wound which include necrotic tissues.

Early diagnosis of non-healing wounds:^{13,14}

The wounds which are not reduced to 50 percent during the first four weeks of therapy and, in fact, have 25 percent reduction after the fourth week in the area of wound (95% CI, 15-35) . In contrast, the healing wounds are healed up to 82 percent during the first four weeks of therapy with the standard therapies (95% CI, 70-94); $p < 0.001$).

Management of the wound:

For the healing and flesh tissue wounds with sizes lower than 3 cm the usual diabetic wound cares including the dressing of the patients with a silver dressing, cleansing of the wounds normal saline solution with

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ROOYAGEN

Autologous PRP-Gel Preparation

This kit should be used by the trained physician

Indication for use:

After wounding in the body, regenerative procedures are activated due to wound healing. In some cases, this repair is not restored in the expected time (less than three months), which is called chronic wounds, and it may take years to improve these ulcers. In some cases, it is never healed, and the progression of the wound leads to amputation of the tissue and ^{1and3}

There are several reasons for chronic wounds most of them are vascular (artery and vein) wounds, diabetic wounds, and pressure wounds resulting in tissue ischemia sequentially. Of these wounds, the healing of diabetic wounds, due to the intrinsic complexity of the disease and its resulting wounds, are of prime importance. ^{2and4-5}

In spite of the great advancement in diagnosis and healing of diabetes, the chronic ulcers of the foot (diabetic foot) in diabetic patients are still the prominent factor of failing in the healing of these patients. A great number of these patients are stricken by different phases of the foot chronic wounds that should receive pharmaceutical treatment for long periods of time. Unfortunately, due to the lack of successful pharmaceutical treatment later, the surgical treatment will be considered. In addition, in some cases, the removal of the necrotic tissues leads to the loss of a large part of the metatarsus, and in case of advancement, the amputation of the related organ, In a wide range study performed by responsible authorities in Japan and the US on Diabetes Mellitus (DM) patients, it was reported that the numbers of the amputations of secondary lower organs annually were 10000 and 66000 for both of countries, respectively ⁶⁻⁷. Although a precise statistics of the subject has not been published formally in our country, Iran, the evidences show that the situation in Iran is comparable to the global statistics.

Platelets Rich Plasma Gel (PRP gel), as an advanced therapy, is used in repairing both chronic and acute wounds. In the report by the US center