ROOYAGEN

Stem Cell Motivation Kit

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
Reference:

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1. Kenneth Kaushansky, Marshall Lichtman, et al , "Williams Hematology text book, 8th Edition" 2010, part 12, chapter 113,114
2. Ronald Hoffman , Bruce Furie, Philip McGlave, et al Hematology text book: Basic Principles and Practice. 6th edition 2013, chapter 12, chapter 126,127
3. Eugenia Wang. Senescent Human Fibroblasts Resist Programmed Cell Death, and Failure to suppress bcll Is Involved1. CANCER RESEARCH, 1995, 55, 2284-2292
4. G P Dimri, X Lee, G Basile, M Acosta, & et al A biomarker that identifies senescent human cells in culture and in aging skin in vivo. PNAS of USA, vol. 92 no. 20, 9363-9367
5. Park HB, Yang JH, Chung KH. Characterization of the cytokine profile of platelet rich plasma (PRP) and PRP-induced cell proliferation and migration: Upregulation of matrix metalloproteinase-1 and -9 in HaCaT cells. Korean J Hematol. 2011;46:265-273.
6. Food & Drug Administration of USA (www.FDA.gov), Cite: CFR21CFR640.25
7. Hillier Text book of Transfusion medicine Second edition : 2007 : chapter 23 : 308-340
8. International cellular medicine society Guideline for PRP. 2011 Page 3
9. Louis Jacques, MD. Tamara Syrek Jensen, JD. James Rollins, MD, et al. Decision Memo for Autologous Blood-Derived Products for Chronic Non-Healing Wounds (CAG-00190R3) centers for Medicare & Medicaid services in USA, 2012.
10. TREATMENT OF OSTEOARTHRITIS OF THE KNEE evidence based guideline, 2th edition. Adopted by the American Academy of Orthopaedic Surgeons Board of Directors. may, 18, 2013
11. enhanced Effect of Platelet-Rich Plasma Containing a New Carrier on Hair Growth, Dermatol Surg 11;37:1721- 1729 Tokyo, Japan. 2011 published by the American Society for Dermatologic Surgery
12. Isabel Cristina, Antonella Tosti, Professor of Clinical Dermatology, Miami, USA. New investigational drugs for AGA. Expert Opin. Investig. Drugs journal 2013 .
13. Mishra A, Pavelko T (2006). "Treatment of chronic elbow tendinosis with buffered platelet-rich plasma". The American Journal of Sports Medicine 34 (11): 1774-8.
14. Spakova T, Rosocha J, Lacko M, Harvanov D, Gharaibeh A. Treatment of knee joint osteoarthritis with autologous platelet-rich plasma in comparison with hyaluronic acid. Am J Phys Med Rehabil. 2012 May; 91(5):411-7.
15. Marx RE . "Platelet-rich plasma: evidence to support its use". Journal of Oral and Maxillofacial Surgery (2004) 62 (4): 489-96.
16. Alessio Redaelli , et al , Face and neck revitalization with PRP: clinical outcome in a series of 23 consecutively treated patients.. J Drugs Dermatol. 2010 May; 9 (5) : 466-72
17. Li ZJ, Choi HI, Choi DK, Sohn KC, Im M, Seo YJ, et al. Autologous platelet-rich plasma: A potential therapeutic tool for promoting hair growth. Dermatol Surg 2012;38:1040-6.
18. Weibrich G, Kleis WK, Hafner G. Growth factor levels in the platelet-rich plasma produced by 2 different methods: Curasan-type PRP kit versus PCCS PRP system. Int J Oral Maxillofac Implants 2002;17:184-90.

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■ When pulling the upper plasma in the first stage, because of the negative pressure inside the tube, which leads to the plasma return from the syringe to the tube, insert a sterile needle beside the IV cannula into the tube.

- Disposable gloves
- Rack
- Phlebotomy requirements: alcohol 70%, tourniquet,

Storage condition & expire date:

- The Expire date labeled on the Kit unless opened or damaged.
- Optimum condition is Room temperature (22-26°C), should not store less than 8°C or more than 30°C.

Procedure:

RooyaGen will recover a small amount of blood (about **35 CC**) from your arm using a small needle and a specially designed syringe. PRP then applied to the treatment area. The entire PRP production process is usually done in less than 30 minutes.

PRP was prepared using a double-spin method:

- 35ML of blood obtained from healthy donor who had provided informed consent was drawing to syringe containing 5 ML anticoagulant.
- This 40 ML (blood+ anticoagulant) transferred to 4 tubes **B** by adaptor
- These tubes were centrifuged by RooyaGen machine at 1600 RPM* for 10 minutes. (If you have another centrifuge, please see the below table).
- Subsequently, the yellow plasma (containing the platelets and leukocytes) was aspirated using a special syringe then transferred to second stage tubes (2 tubes **C**) and then centrifuged at 3500 RPM** for 6 minutes, yielding a platelet-rich pellet and platelet-poor plasma (PPP).
- The platelet-rich pellet, depended on indication, was re-suspended in 2 ML of plasma (for hair), 3 ML of plasma (for orthopedic), 4 ML of plasma (for skin) and the resulting suspension used as PRP.
- According to site, PRP inject with 27 g or 30 g needle injection into the true cortex, really sub-cortical or periosteum plane before.
- For musculoskeletal 21G to 25G needle injection in target tissue.

Note: * and ** due to variation in centrifuge, please set up your centrifuge according to below table:

Brand	Model	Status	Country	1ST(RPM)	2ST(PRM)
Rooyagen	-	Swing	Iran	1600	3500
Hettich	EBA 20	Fix	Germany	1200	3300
Hettich	EBA 270	Swing	Germany	2700	4000
Hettich	Rotofix 32	Swing	Germany	1600	3500
Pars azma	-	Swing	Iran	1500	3400
Jenus	TDL-4	Fix	China	1000	2700
Sehand	-	Swing	Iran	1600	3500
Labtron	-	Swing	Iran	1600	3500
Selecta La	-	Fix	S. Korea	900	2400
Selecta La	-	Swing	S.Korea	1300	3300

Side effect & cautions:

PRP uses your body's own natural properties to treat your injury. Side effects utilizing PRP systems are very uncommon unless usage of IVD class device and without authorized Kit which made of non-medical plastic (industrial material is not suitable for PRP kit). If contraindication was seen don't saw any side effects But slightly bruising occurs in some case (special in eye lids) that could improvement in 3 to 5 days.

Injection protocol & results:

About optimum dose for non-transfusion PRP therapy despite the controversy between publications, the mean of acceptable articles have a partial agreement about 3.5 time rich more than baseline for skin^{11,12}, 5 time rich more than base line for musculoskeletal^{13,14} and 7 time rich more than baseline for hair growth^{15,16} which should be confirmed by laboratory cell analyzer (CBC test).

For skin: 2 to 3 injection with 2 month interval (between 3 to 6 weeks) and one annually remember injection was needed. Definite results was seen 4 month after first injection, probably saw the results 3 to 4 weeks after that but the maximum saw in 4 month.

For hair: 2 to 4 injection with 3 to 6 week interval and one annually remember injection was needed. Definite results were seen 5 month after first injection, probably saw the results such as hair-loose disconnect 2 to 3 weeks after that but the maximum saw in 4 to 5 month.

For tendon and ligaments: 2 injection with 3 to 6 weeks interval and one annually remember injection was needed. Definite results was seen during the first month(probably saw the results 7 to 10 day after that but the maximum saw in 3 weeks).

Indication:

New research and technology has expanded the application of this therapy for use in orthopedics, cosmetic, plastic surgery, dentist, wound healing and...⁸. At now several PRP Kit was approved by FDA (Food & Drug administration in USA) and CE (Council Europa) and this therapy was approved by CMS⁹ (centers for Medicare & Medicaid services in USA). Recently in five-year report by AAOS (American academy of Orthopedic surgery) PRP approved as one of the non-invasive methods that more effective than Hyaluronic Acid & steroids.¹⁰ More than 5200 published clinical studies can be found on NCBI⁸, investigating the effectiveness of Autologous Platelet Rich Plasma in various medical disciplines such as **musculoskeletal**, **sports medicine**, **wound healing**, **surgery** and **aesthetics**, but definitely must usage of PRP including:

Fine line and wrinkle	Textural improvement
Dull dry skin	Hair growth stimulation (for thinning hair)
Full face rejuvenation	Eye area (both upper & lower eyelids)
Décolletage	Back of hands
Diabetic & pressure wound	Long standing problems with tendonitis & arthritis
Plastic surgery	And ...

Contraindication:

According to ICMS guideline- 2011 (International Cellular Medicine Society)

ABSOLUTE CONTRAINDICATIONS:	RELATIVE CONTRAINDICATIONS:
<ul style="list-style-type: none">• Platelet dysfunction syndrome• Critical thrombocytopenia (<50000)• Hemodynamic instability• Septicemia• Local infection at the site of the procedure• Patient unwilling to accept risks• Do not use in pregnant women and children	<ul style="list-style-type: none">• Consistent use of NSAIDs within 48 hours of procedure• Corticosteroid injection at treatment site within 1 month• Systemic use of corticosteroids within 2 weeks• Tobacco use• Recent fever or illness• Cancer- especially hematopoietic or of bone• HGB < 10 g/dl• Platelet count < 10⁵/ul

Cautions:

- Blood drawing without any stress on vein in limited time (about 3-5 minute) by trained person
- Gently act in time period for Platelet separation because rigid working could stimulate of platelet and due to growth factor releasing or apoptosis initiation.
- To be sure from centrifuge balanced and sharp border between blood and plasma, if don't see loose the centrifuge table, tube unbalanced and finally centrifuge abrupt to be cautioned.
- Approximately four to seven times the mean levels is generally regarded as the Therapeutically effective concentration of PRP.^{17,18} and this increased richness should be confirmed by CBC test.
- For achievement to good results should use the rich in any cases, for example about **3.5 time rich for skin**^{11,12}, about **5 time rich for musculoskeletal**^{13,14} and **7 time rich for hair growth**^{15,16} that confirm by laboratory cell analyzer. (CBC test)
- When transfer the plasma from stage 1 to stage 2, exactly relax and prevention from RBCs transfer, (the acceptance level of RBCs in PRP should be less than 2.0%).
- Needle gauge for skin & hair 27G to 30G
- Needle gauge for tendon and ligament 23G to 25G

Kit contents:

STAGE I

Butterfly Needle (21g) 1X
Syringe For Blood Collection (50ml) 1X
Anti Coagulant and Separator (6.5ml) 1X
Platelet Non-Adhesive Tube (10ml) 4X
Vacutainer Adaptor For Transfer 1X

STAGE II

Needle For Transfer (16g) 2X
PRP Preservative Tube (10ml) 2X
Syringe For Supernatant Transfer (5ml) 2X
Needle For Transfer (21g) 2X
Insulin Syringe (1ml) 1X
Needle For Injection (27-30g) For Cosmetic 4X


Equipment is not available in Kit:

- Centrifuge for 16*100 tube (swing was preferred)
- Safety Box

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For autologous plasma rich platelet (PRP) separation

Introduction:

Blood including two major content: the first is cellular content which including RBC, WBC and Platelet and second is plasma which contain water (about 93%) , various proteins (Albumin, anti-body , coagulation factors,...) trace elements (Iron , Zinc, calcium,... , vitamins,...).¹

Platelet is a part of cellular component which has 10 day half time (9-11 day) and it was about 6% of blood count and/or 0.4% of blood volume. PRP is a concentration of platelets and growth factors created from a small amount of your own blood. It has two function, the first is cofactor for blood coagulation to prohibition of bleeding and second Increased levels of growth factors improve signaling and recruitment of cells to an injury site and optimize the environment for healing.²

Healing after an injury involves a well-orchestrated and complex series of events where proteins in the blood act as messengers to regulate the entire process. Many proteins involved in the healing process are derived from small cell fragments in the blood called platelets. Same as injury , in the middle decade of live the stem cell going to senescence phase and don't response to normal level of factors , however they needed to enrichment of growth factors more than baseline level. The platelets become activated and start together at the injury site to release beneficial proteins called growth factors. This is the beginning of the rejuvenation process.^{3,4,5}

The major complication in PRP therapy is various separation techniques which shall be cautioned the FDA⁶ requirements by Cite No: 21CFR640.25, in this letter in addition swirling, that should be see, the Platelet count should be analyze by CBC test in laboratory at least for 1% of patient. Finally, only two-round centrifuge method documented in Hematology & Transfusion medicine text book for PRP separation.⁷